



ELK Biotechnology
For research use only.

High Density Lipoprotein – Cholesterol (HDL-C) Assay Kit

Microplate Method




CAT/NO.: BC029 96T

1. Reagents Composition and Preparation

Reagents	Format	Composition	Concentration	Storage Condition
Reagent I	1 Bottle×18mL	Good's Buffer	50mmol/L	2-8°C Avoid Light
		N-Ethyl-N-(2-Hydroxy-3-Sulfopropyl)-3-Methylaniline Sodium Salt (Toos)	1mmol/L	
		Magnesium Chloride Hexahydrate	15mmol/L	
		Cholesterol Oxidase (CO)	≥3kU/L	
		Peroxidase	≥5kU/L	
Reagent II	1 Bottle×6mL	Good's Buffer	50mmol/L	
		4-Amino Antipyrine	0.2mmol/L	
		Magnesium Chloride Hexahydrate	15mmol/L	
		Cholesteryl Esterase (CE)	≥3kU/L	
		Surfactants	0.1%	
Standard	1 Bottle	Cholesterol	See label	

2. Principle of Measurement

1. First Reaction



Note 1. Δ^4 represents Cholest-4-en-3-one



ELK Biotechnology

For research use only.

3. Procedures of Measurement

I. Sample Pre-Treatment

1. Serum (Plasma): Can be used for further measurement directly. The serum should be diluted with physiological saline if the concentration of LDL-C in serum exceeds the concentration range allowed.
2. Growth Medium Sample: Sample medium should be extracted and separated by centrifugation at 2,500rpm for 10 minutes. Supernatant is desired.
3. Tissue Sample: Samples should be weighed precisely and for every 1g of the sample, 9mL homogenate medium is mixed with sample in ice water bath. The homogenate treated with centrifugation at 2500rpm for 10 minutes would be separated and supernatant is extracted for further use.

Note 1. The homogenate medium is 0.1mmol/L phosphate buffer (pH=7.4) or physiological saline for tissue samples without high fat part.

Note 2. For high fat or partial high fat sample, the homogenate medium used is absolute alcohol.

4. Sample Cells

i. Cell Harvesting

Cell suspension is extracted and centrifuged with 1000rpm for 10 minutes. Supernatant should be discarded and isotonic buffer (0.1mol/L, pH 7-7.4 phosphate buffer recommended) is used to wash the precipitates for once or twice. Suspension is separated by centrifugation for 10 minutes and precipitates should be kept for further treatment.

ii. Cell Disruption

0.2 to 0.3 mL homogenate medium should be added and homogenized (0.1mol/L, pH 7-7.4 phosphate buffer or physiological saline recommended). Cells can be disrupted with ultrasonic disruptor in ice water bath. The ultrasonic disruptor is 300W and the disruption is repeated for 3-5 times with the period of 3-5 seconds and interval of 30 seconds. The homogenate can be measured directly. Also, Lysate can be used directly for the later measurement (Triton X-100, 1-2% recommended for the splitting period of 30-40 minutes).

Note that cell density in sample is recommended to be no less than 1 million/ml. Microscope can be used to observe the completeness of cell disruption.



ELK Biotechnology
For research use only.

4. Measurement Procedures

Compositions (μ l)	Blank	Standard	Sample
Distilled Water	2.5		
1.8mM Standard Solution		2.5	
Sample			2.5
Reagent I	180	180	180
Mix thoroughly and warm at 37° C for 5 min. Use enzyme-labeled instrument to measure absorbances A ₁ at 600nm.			
Reagent II	60	60	60

Mix thoroughly and warm at 37° C for 5-10 min. Use enzyme-labeled instrument to measure absorbances A₂ at 600nm.

5. Calculation Formula and Example

1. Serum or other Liquid Sample

$$HDL - C = \frac{(A_{2,Sample} - A_{1,Sample}) - (A_{2,Blank} - A_{1,Blank})}{(A_{2,Standard} - A_{1,Standard}) - (A_{2,Blank} - A_{1,Blank})} \times \frac{C_{Standard}}{1.8mM}$$

Example: 2.5μL human serum was measured and absorbance values are listed below:

$A_{1,Blank} = 0.000$, $A_{2,Blank} = 0.010$, $A_{1,Standard} = 0.048$, $A_{2,Standard} = 0.139$, $A_{1,Sample} = 0.042$ and $A_{2,Sample} = 0.066$. The standard reagent's concentration is 1.8mM.

The HDL-C concentration is

$$HDL - C = \frac{(A_{2,Sample} - A_{1,Sample}) - (A_{2,Blank} - A_{1,Blank})}{(A_{2,Standard} - A_{1,Standard}) - (A_{2,Blank} - A_{1,Blank})} \times \frac{C_{Standard}}{1.8mM}$$

$$= \frac{(0.066 - 0.042) - (0.010 - 0.000)}{(0.139 - 0.048) - (0.010 - 0.000)} \times 1.8 = 1.273mM$$

Example: 2.5μL serum taken from rat was measured and absorbance values are listed below:

$A_{1,Blank} = 0.010$, $A_{2,Blank} = 0.000$, $A_{1,Standard} = 0.048$, $A_{2,Standard} = 0.139$, $A_{1,Sample} = 0.027$ and $A_{2,Sample} = 0.005$. The standard reagent's concentration is 1.8mM.

$$HDL - C = \frac{(A_{2,Sample} - A_{1,Sample}) - (A_{2,Blank} - A_{1,Blank})}{(A_{2,Standard} - A_{1,Standard}) - (A_{2,Blank} - A_{1,Blank})} \times \frac{C_{Standard}}{1.8mM}$$

$$= \frac{(0.005 - 0.027) - (0.000 - 0.010)}{(0.139 - 0.048) - (0.010 - 0.000)} \times 1.8 = 0.2744mM$$

2. Tissue or Cell Sample

For Microplate Reader:



ELK Biotechnology

For research use only.

$$HDL - C = \frac{(A_{2,Sample} - A_{1,Sample}) - (A_{2,Blank} - A_{1,Blank})}{(A_{2,Standard} - A_{1,Standard}) - (A_{2,Blank} - A_{1,Blank})} \times \frac{C_{Standard}}{1.8mM} \div C(mgprot/g)$$

Example: 2.5 μ L homogenate from mouse's liver was measured and absorbance values are listed below:

$A_{1,Blank} = 0.000$, $A_{2,Blank} = 0.010$, $A_{1,Standard} = 0.048$, $A_{2,Standard} = 0.139$, $A_{1,Sample} = 0.066$ and $A_{2,Sample} = 0.040$. The standard reagent's concentration is 1.8mmol/L. Also, the protein concentration is measured to be 12.0121 g/L.

$$\begin{aligned} HDL - C &= \frac{(A_{2,Sample} - A_{1,Sample}) - (A_{2,Blank} - A_{1,Blank})}{(A_{2,Standard} - A_{1,Standard}) - (A_{2,Blank} - A_{1,Blank})} \times C_{Standard} \div \rho(mgprot/mg) \\ &= \frac{(0.066 - 0.060) - (0.010 - 0.000)}{(0.139 - 0.048) - (0.010 - 0.000)} \times 1.8 \div 12.0121 = 0.0292mmol/g \end{aligned}$$

6.Performance Index

1. Blank tube absorbance value \leq 0.005.
2. Sensitivity: The absorbance value difference ΔA is 0.087-0.153 for the 1.3mM standard solution.
3. Restriction: Concentration range permitted is 0.65-3.8mmol/L with $r^2 > 0.995$.
4. Accuracy: The coefficient of variation \leq 3% and relative deviation among batches \leq 5%
5. Stability: Reagent in the original package should be preserved at 2-8 $^{\circ}$ C without light-struck and is stable within 12 months. Unsealed packages can be preserved for one month under same condition.

7.Notes

1. This assay kit is designed strictly for scientific research. Usage for clinical trials or oral taken is prohibited.
2. Please dilute the sample with physiological saline when the concentration exceeds the range permitted.
3. Reagents can be contaminated by glucose or cholesterol assay reagent.
4. The amount of reagent and sample tested can be adjusted according to the instructions of the automatic biochemical analyzer.