

## **ELK Biotechnology**

GFAP Mouse mAb Catalog NO.: EM1059 For research use only.

Overview


Product name	GFAP Mouse Monoclonal antibody
Source	Mouse
Applications	WB IHC
Species reactivity	Rat Mouse
Recommended dilutions	WesternBlot:1/2000-5000 Immunohistochemistry:1/200-500 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic Peptide
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
lsotype	lgG1
Clonality	Monoclonal
Concentration	1 mg/ml
Observed band	45kDa
GenelD (Human)	2670
Human Swiss-Prot No.	P14136
Cellular localization	Cytoplasm. Associated with intermediate filaments
Alternative Names	Glial Fibrillary Acidic Protein ;Intermediate Filament Protein
Background	GFAP is a member of the class III intermediate filament protein family. It is heavily and specifically expressed in astrocytes and certain other astroglia in the central nervous system in satellite cells in peripheral ganglia and in non myelinating Schwann cells in peripheral nerves. In addition neural stem cells frequently strongly express GFAP. Antibodies to GFAP are therefore very useful as markers of astrocytic cells. In addition many types of brain tumor presumably derived from astrocytic cells heavily express GFAP.
https://www.elkbiotech.com	T:86-27-59760950 E: elkbio@elkbiotech.com

GFAP is also found in the lens epithelium Kupffer cells of the liver in some cells in salivary tumors and has been reported in erythrocytes.

Western blot analysis of Rat Brain Tissue with GFAP mAb diluted at:5000.

Western blot analysis of Mouse Brain Tissue with GFAP mAb diluted at:2000.

Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using GFAP (EM1059) Mouse mAb diluted at:500.



94KD 66KD