

ErbB-2/HER-2 Mouse mAb Catalog NO.: EM1079 For research use only.

## Overview

**Product name** ErbB-2/HER-2 Mouse Monoclonal antibody

**Source** Mouse

Applications WB IHC

Species reactivity Human

Recommended dilutions WesternBlot:1/1000

Immunohistochemistry:1/200-500

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Synthetic Peptide

**Species** Human

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

lsotype lgG1

**Clonality** Monoclonal

Concentration 1 mg/ml

Observed band 138kDa

GenelD (Human) 2064

Human Swiss-Prot No. P04626

**Cellular localization** Cytoplasm. Nucleus and Cell membrane

Alternative Names CD340 HER 2/neu HER2 MLN9 MLN19 NEU NGL p185 p185erbB2 Proto

oncogene c ErbB 2 Proto oncogene Neu TKR1


Background ErbB 2 is a receptor tyrosine kinase of the ErbB 2 family. It is closely related

instructure to the epidermal growth factor receptor. ErbB 2 oncoprotein is detectable in a proportion of breast and other adenocarconomas as well as transitional cell carcinomas. In the case of breast cancer expression


determined by immunohistochemistry has been shown to be associated with

poor prognosis.

https://www.elkbiotech.com T:86-27-59760950 E: elkbio@elkbiotech.com



Western blot analysis of) Hela 2) 293 3) Jurkatwith ErbB-2 Mouse mAb diluted at:2000.



Immunohistochemical analysis of paraffin-embedded Human Breast Caricnoma using ErbB-2/HER-2 (EM1079) Mouse mAb diluted at:500.