

ELK Biotechnology

alpha Lactalbumin Mouse mAb Catalog NO.: EM1084

For research use only.

Overview

Product name alpha Lactalbumin Mouse Monoclonal antibody

Source Mouse

Applications WB

Species reactivity Human Bovine

Recommended dilutions WesternBlot:1/1000

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Synthetic Peptide

Species Human

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

Isotype IgG1

Clonality Monoclonal

Concentration 1 mg/ml

Observed band 16kDa

GenelD (Human) 3906


Human Swiss-Prot No. P00709

Cellular localization Secreted

Alternative Names N/A

Background In primates alpha-lactalbumin expression is upregulated in response to the

hormone prolactin and increases the production of lactose. α -Lactalbumin forms the regulatory subunit of the lactose synthase (LS) heterodimer and β -14-galactosyltransferase (beta4Gal-T1) forms the catalytic component.



Western blot analysis of) Human Milk 2) Milk with alpha Lactalbumin Mouse mAb diluted at: 3000.