

ELK Biotechnology

Rabbit IgG(H+L) Mouse mAb

Catalog NO.: EM1167 For research use only.

Overview

Product name Rabbit IgG(H+L) Mouse Monoclonal antibody

Source Mouse

Applications WB IHC

Species reactivity Rat

Recommended dilutions WesternBlot:1/5000-20000

Immunohistochemistry:1/200-500

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Recombinant Protein

Species Rabbit

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

lsotype lgG1

Clonality Monoclonal

Concentration 1 mg/ml

Observed band 2555kDa

GenelD (Human) N/A


Human Swiss-Prot No. N/A

Cellular localization N/A

Alternative Names N/A

Background The antibody reacts with whole molecule rabbit IgG. The antibody may

cross-react with immunoglobulins from other species.



Western blot analysis of Reduced Rabbit IgG H+L with EM1167 diluted at: 10000.