

RFP-Tag Mouse mAb Catalog NO.: EM1170 For research use only.

## Overview

Product name RFP-Tag Mouse Monoclonal antibody

**Source** Mouse

Applications WB IHC IP

Species reactivity N/A

Recommended dilutions WesternBlot:1/5000

Immunohistochemistry:1/200 Immunoprecipitation:1/200

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Recombinant Protein

Species N/A

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

lsotype lgG1

**Clonality** Monoclonal

Concentration 1 mg/ml

Observed band N/A

GenelD (Human) N/A

Human Swiss-Prot No. N/A

Cellular localization N/A

Alternative Names DsRed Red fluorescent protein Red fluorescent protein drFP583

**Background** The red fluorescent protein cloned from Discosoma coral (DsRed or

drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companionor substitute for the green fluorescent

protein(GFP) from the Aequorea jellyfish.

https://www.elkbiotech.com

T:86-27-59760950

E: elkbio@elkbiotech.com



Western blot analysis of recombinant RFP protein with EM1170 diluted at: 10000.