

PTEN (2C10) Mouse mAb Catalog NO.: EM1308 For research use only.

Overview

Product name PTEN (2C10) Mouse Monoclonal antibody

Source Mouse

Applications IHC

Species reactivity Human

Recommended dilutions Immunohistochemistry:1/100-200

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Synthetic Peptide

Species Human

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

Isotype IgG1

Clonality Monoclonal

Concentration 1 mg/ml

Observed band 54kDa

GenelD (Human) 5728

Human Swiss-Prot No. P60484


Cellular localization Cytoplasm Nucleus Secreted

Alternative Names BZSMHAMMMAC1TEP1phosphatase and tensin homolog

Background PTEN (phosphatase and tensin homologue deleted on chromosome ten)

also referred to as MMAC (mutatedIn multiple advanced cancers) phosphatasels a tumor suppressorImplicatedIn a wide variety of human

cancers.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using PTEN (EM1308) Mouse mAb diluted at:200.