


ELK Biotechnology ERK1 (5E9) Mouse mAb Catalog NO.: EM1323 For research use only.

Overview

Product name	ERK1 (5E9) Mouse Monoclonal antibody
Source	Mouse
Applications	IHC
Species reactivity	Human Rat Mouse
Recommended dilutions	Immunohistochemistry:1/100-200 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1 mg/ml
Observed band	44kDa
GenelD (Human)	5594
Human Swiss-Prot No.	P27361
Cellular localization	Nucleus
Alternative Names	MAPK Mitogen-activated protein Kinase PRKM2MK01
Background	Mitogen-activated protein kinases (MAPKs) are a widely conserved family of serine/threonine protein kinasesInvolvedIn many cellular programs such as cell proliferation differentiation motility and death. The p44/42 MAPK (Erk1/2) signaling pathway can be activatedIn response to a diverse range of extracellular stimuliIncluding mitogens growth factors and cytokines andIs anImportant targetIn the diagnosis and treatment of cancer.



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using ERK1 (EM1323) Mouse mAb diluted at:200.

Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using ERK1 (EM1323) Mouse mAb diluted at:200.