

MLKL (8H7) Mouse mAb Catalog NO.: EM1328 For research use only.

Overview

Product name MLKL (8H7) Mouse Monoclonal antibody

Source Mouse

Applications IHC

Species reactivity Human

Recommended dilutions Immunohistochemistry:1/100-200

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Synthetic Peptide

Species Human

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

Isotype IgG1

Clonality Monoclonal

Concentration 1 mg/ml

Observed band 54kDa

GenelD (Human) 197259

Human Swiss-Prot No. Q8NB16


Cellular localization Cell membrane Cytoplasm Membrane

Alternative Names Mixed lineage kinase domain like protein FLJ34389


Background Mixed lineage kinase domain-like protein (MLKL)Is a pseudokinase that

wasIdentified as downstream target of RIP3In the necroptosis pathway. During necroptosis RIP3Is phosphorylated at Ser227 which recruits MLKL

and leads tolts phosphorylation at Thr357 and Ser358.



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using MLKL (EM1328) Mouse mAb diluted at:200..



Immunohistochemical analysis of paraffin-embedded Human Skin Tissue using MLKL (EM1328) Mouse mAb diluted at:200.