

ELK Biotechnology

HA-Tag Rabbit pAb Catalog NO.: EA002 For research use only.

Overview

Product name	HA-Tag Rabbit polyclonal antibody
Source	Rabbit
Applications	WB
Species reactivity	N/A
Recommended dilutions	WesternBlot:1/5000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic Peptide
Species	N/A
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
lsotype	lgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	N/A
GenelD (Human)	N/A
Human Swiss-Prot No.	N/A
Cellular localization	N/A
Alternative Names	N/A
Background	Human influenza hemagglutinin (HA) is a surface glycoprotein required for the infectivity of the human virus. The HA tag is derived from the HA- molecule corresponding to amino acids 98-106 has been extensively used as a general epitope tag in expression vectors. Many recombinant proteins have been engineered to express the HA tag, which does not appear to interfere with the bioactivity or the biodistribution of the recombinant protein.HA tag antibody is a highly sensitive and affinity monoclonal antibody applicable to HA-tagged fusion protein detection.
https://www.elkbiotech.com	T:86-27-59760950 E: elkbio@elkbiotech.com 3

