

ELK Biotechnology

MeCP2 Rabbit pAb Catalog NO.: EA020 For research use only.

Overview

Product name	MeCP2 Rabbit polyclonal antibody
Source	Rabbit
Applications	WB
Species reactivity	Human
Recommended dilutions	WesternBlot:1/2000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic Peptide
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
lsotype	lgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	53kDa
GenelD (Human)	4204
Human Swiss-Prot No.	P51608
Cellular localization	Nucleus
Alternative Names	Methyl CpG Binding Protein 2 ;RTS; AUTSX3; MRX16; MRX79; PPMX; RTT; Rett Syndrome; Mental Retardation,X-linked 16
Background	Methyl-CpG-binding protein 2(MeCP2), whose exact function is unclear. It appears to help regulate gene activity (expression) by modifying chromatin, the complex of DNA and protein that packages DNA into chromosomes. The MeCP2 protein usually regulates genes involved in brain function, even though this protein is found throughout the body.Within the brain, the MeCP2 protein is important for the function of nerve cells (neurons) and is present in high levels in mature neurons. This protein likely plays a role in
https://www.elkbiotech.com	T:86-27-59760950 E: elkbio@elkbiotech.com 3(

maintaining connections (synapses) between neurons, where cell-to-cell communication occurs. Many of the genes that are known to be regulated by the MeCP2 protein play a role in normal brain function, particularly the maintenance of synapses.

Western blot analysis of Hela with MeCP2 Rabbit pAb diluted at 1, 1:2,000 2, 1:5,000

