



ELK Biotechnology


EFHD1 Rabbit pAb

Catalog NO.: EA021

For research use only.

Overview

Product name	EFHD1 Rabbit polyclonal antibody
Source	Rabbit
Applications	WB
Species reactivity	Human
Recommended dilutions	WesternBlot:1/1000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	27kDa
GeneID (Human)	80303
Human Swiss-Prot No.	Q9BUP0
Cellular localization	Membrane, Mitochondrion, Mitochondrion inner membrane
Alternative Names	MST133, MSTP133, PP3051, Swiprosin 2, SWS2
Background	EFHD1 is an EF-hand domain-containing protein that displays increased expression during neuronal differentiation.



Western blot analysis of 293T cell lysate with EFHD1 rabbit polyclonal antibody diluted at 1:2,000.