

GAPDH Rabbit pAb Catalog NO.: EA027 For research use only.

Overview

Product name GAPDH Rabbit polyclonal antibody

Source Rabbit

Applications WB, IHC

Species reactivity Zebrafish

Recommended dilutions WesternBlot:1/5000-10000

Immunohistochemistry:1/200 Immunofluorescence:1/200

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Recombinant Protein

Species Human

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

Isotype IgG

Clonality Polyclonal

Concentration 1 mg/ml

Observed band 37kDa

GenelD (Human) 2597

Human Swiss-Prot No. P04406

Cellular localization N/A

Alternative Names G3PD, G3PDH, GAPD, Glyceraldehyde 3 phosphoate dehycdrogenase,

MGC102544

Background Glyceraldehyde 3 phosphate dehydrogenase (GAPDH) is well known as

one of the key enzymes involved in glycolysis. GAPDH is constitutively expressed in almost all tissues at high levels, therefore antibodies against GAPDH are useful as loading controls for Western Blotting. Some physiological factors, such as hypoxia and diabetes, increase GAPDH

https://www.elkbiotech.com T:86-27-59760950

E: elkbio@elkbiotech.com

expression in certain cell types.



