

## ECFP-Tag Rabbit pAb Catalog NO.: EA028

For research use only.

## Overview

Product name ECFP-Tag Rabbit polyclonal antibody

**Source** Rabbit

Applications WB

Species reactivity N/A

Recommended dilutions WesternBlot:1/2000-5000

Immunoprecipitation:1/200

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Recombinant Protein

Species N/A

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

**Isotype** IgG

**Clonality** Polyclonal

Concentration 1 mg/ml

Observed band N/A

GenelD (Human) N/A

Human Swiss-Prot No. N/A

Cellular localization N/A

Alternative Names Cyan fluorescent protein, Enhanced Cyan fluorescent protein

**Background**Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent

protein (GFP) originally derived from the jellyfish Aequorea victoria.



Western blot analysis of ECFP recombinant protein with EA028 diluted at 1:5,000.