

ELK Biotechnology

mStrawberry-Tag Rabbit pAb

Catalog NO.: EA034 For research use only.

Overview

Product name mStrawberry-Tag Rabbit polyclonal antibody

Source Rabbit

Applications WB

Species reactivity N/A

Recommended dilutions WesternBlot:1/2000-5000

Immunohistochemistry:1/200 Immunofluorescence:1/200 Immunoprecipitation:1/200

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Recombinant Protein

Species N/A

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

Isotype IgG

Clonality Polyclonal

Concentration 1 mg/ml

Observed band N/A

GenelD (Human) N/A

Human Swiss-Prot No. N/A

Cellular localization N/A

Alternative Names mStrawberry

Background mStrawberry is a bright, monomeric red fluorescent protein which was

developed in Dr. Roger Tsien's lab by directed mutagenesis of mRFP.



Western blot analysis of mStawberry recombinant protein with EA034 diluted at 1:5,000.