



ELK Biotechnology


TLR13 Rabbit pAb

Catalog NO.: EA182

For research use only.

Overview

Product name	TLR13 Rabbit polyclonal antibody
Source	Rabbit
Applications	IHC
Species reactivity	Human, Rat, Mouse
Recommended dilutions	Immunohistochemistry:1/200-500 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	N/A
GeneID (Human)	279572
Human Swiss-Prot No.	Q6R5N8
Cellular localization	Plasma membrane
Alternative Names	Toll like receptor 13 antibody
Background	TLR13 is found in mouse and participates in the innate immune response to microbial agents. It acts via MYD88 and TRAF6, leading to NF-kappa-B activation, cytokine secretion and the inflammatory response.



Immunohistochemical analysis of paraffin-embedded Mouse Spleen Tissue using TLR13 (EA182) Rabbit pAb diluted at 1:500.