



ELK Biotechnology

GAL4 DNA Binding Domain Rabbit pAb

Catalog NO.: EA202

For research use only.

Overview

Product name	GAL4 DNA Binding Domain Rabbit polyclonal antibody
Source	Rabbit
Applications	WB
Species reactivity	N/A
Recommended dilutions	WesternBlot:1/10000-50000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	N/A
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	N/A
GeneID (Human)	N/A
Human Swiss-Prot No.	N/A
Cellular localization	N/A
Alternative Names	N/A
Background	The GAL4 protein of <i>Saccharomyces cerevisiae</i> is one of the most thoroughly characterized transcriptional activators. Since the N-terminal 147 amino acid residues of GAL4 are sufficient to mediate specific and strong binding to DNA.



Western blot analysis of Recombinant BD protein with GAL4 Binding Domain Rabbit pAb EA202 diluted at 1:20,000.