



ELK Biotechnology

GALR2 Rabbit pAb

Catalog NO.: EA291

For research use only.

Overview

Product name	GALR2 Rabbit polyclonal antibody
Source	Rabbit
Applications	WB
Species reactivity	Human, Rat, Mouse
Recommended dilutions	WesternBlot:1/1000-2000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic Peptide
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	42kDa
GeneID (Human)	8811
Human Swiss-Prot No.	O43603
Cellular localization	Cell membrane, Membrane
Alternative Names	Galanin receptor 2,MGC125983
Background	GALR2 Receptor for the hormone galanin and for GALP. The activity of this receptor is mediated by G proteins that activate the phospholipase C/protein kinase C pathway (via G(q)) and that inhibit adenylyl cyclase (via G(i)).

