

SLC12A4 Rabbit pAb Catalog NO.: EA303 For research use only.

## Overview

Product name SLC12A4 Rabbit polyclonal antibody

**Source** Rabbit

Applications WB, IHC

Species reactivity Human, Rat, Mouse

Recommended dilutions WesternBlot:1/1000-2000

Immunohistochemistry:1/100-200

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Synthetic Peptide

Species Human

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

**Isotype** IgG

**Clonality** Polyclonal

Concentration 1 mg/ml

Observed band 150kDa

GenelD (Human) 6560

Human Swiss-Prot No. Q9UP95


Cellular localization Membrane

Alternative Names Electroneutral potassium chloride cotransporter 1,hKCC1,KCC1,


s12A4,FLJ40489

**Background** Mediates electroneutral potassium-chloride cotransport when activated by

cell swelling. May contribute to cell volume homeostasis in single cells.



Western blot analysis of 1) Rat Brain Tissue, 2)Mouse Brain Tissue, 3) Hela with SLC12A4 Rabbit pAb EA303 diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Rat BrainTissue using SLC12A4 (EA303) Rabbit pAb diluted at 1:200.

Immunohistochemical analysis of paraffin-embedded Mouse BrainTissue using SLC12A4 (EA303) Rabbit pAb diluted at 1:200.