



ELK Biotechnology


Kv1.5 Rabbit pAb

Catalog NO.: EA309

For research use only.

Overview

Product name	Kv1.5 Rabbit polyclonal antibody
Source	Rabbit
Applications	WB
Species reactivity	Human
Recommended dilutions	WesternBlot:1/1000-2000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic Peptide
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	67kDa
GeneID (Human)	3741
Human Swiss-Prot No.	P22460
Cellular localization	Membrane
Alternative Names	KCNA5,PCN5,Potassium channel 1, voltage gated potassium channel protein kv1.5
Background	N/A



Western blot analysis of 293T with KV1.5 Rabbit pAb EA309 diluted at 1:2,000.