



ELK Biotechnology

DUSP6 Rabbit pAb

Catalog NO.: EA315

For research use only.

Overview

Product name	DUSP6 Rabbit polyclonal antibody
Source	Rabbit
Applications	WB
Species reactivity	Human, Rat, Mouse
Recommended dilutions	WesternBlot:1/1000-5000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic Peptide
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	42,44kDa
GeneID (Human)	1848
Human Swiss-Prot No.	Q16828
Cellular localization	Cytoplasm
Alternative Names	Dual specificity phosphatase 6, HH19, MKP3
Background	N/A

