

Cleaved-Caspase-3 p12 (D175) rabbit pAb

Cat No.:ES1001

For research use only

Overview

Product Name	Cleaved-Caspase-3 p12 (D175) rabbit pAb Rabbit
Annlications	
Species Cross-Reactivity	Human:Rat:Mouse:
Becommended dilutions	WB 1.500-2000 IHC-n 1.50-300 IE 1.50-300
Immunogen	The antiserum was produced against synthesized
	peptide derived from human Caspase 3. AA
	range:157-206
Specificity	Cleaved-Caspase-3 p12 (D175) Polyclonal Antibody
	detects endogenous levels of fragment of activated
	Caspase-3 p12 protein resulting from cleavage
	adjacent to D175.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$. Avoid repeated freeze-thaw cycles.
Protein Name	Caspase3
Gene Name	CASP3
Cellular localization	Cytoplasm.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	12 35kD
Human Gene ID	836
Human Swiss-Prot Number	P42574
Alternative Names	CASP3; CPP32; Caspase-3; CASP-3; Apopain;
	Cysteine protease CPP32; CPP-32; Protein Yama;
	SREBP cleavage activity 1; SCA-1
Background	This gene encodes a protein which is a member of
	the cysteine-aspartic acid protease (caspase) family.
	Sequential activation of caspases plays a central role
	in the execution-phase of cell apoptosis. Caspases



exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein cleaves and activates caspases 6, 7 and 9, and the protein itself is processed by caspases 8, 9 and 10. It is the predominant caspase involved in the cleavage of amyloid-beta 4A precursor protein, which is associated with neuronal death in Alzheimer's disease. Alternative splicing of this gene results in two transcript variants that encode the same protein. [provided by RefSeq, Jul 2008],