



ELK Biotechnology

PKD1 Rabbit pAb

Catalog Number: ES10021

For research use only.

Overview

Product Name	PKD1 polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human
Recommended dilutions	Immunohistochemistry:1/50-1/300
Immunogen	Synthesized peptide derived from human protein . at AA range; 1370-1450
Specificity	PKD1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	PKD1
Gene Name	PKD1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	473kDa
Human Gene ID	5310
Human Swiss-Prot Number	P98161
Alternative Names	