



ELK Biotechnology

BID Rabbit pAb
Catalog Number: ES1012
For research use only.

Overview

Product Name	BID polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human BID. at AA range; 20-100
Specificity	BID Polyclonal Antibody detects endogenous levels of BID protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	BID
Gene Name	BID
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	21kDa
Human Gene ID	637
Human Swiss-Prot Number	P55957
Alternative Names	BID; BH3-interacting domain death agonist; p22 BID; BID