



# RAB30 rabbit pAb

Cat No.:ES10121

For research use only


## Overview

---

<b>Product Name</b>	RAB30 rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB;ELISA
<b>Species Cross-Reactivity</b>	Human;Mouse
<b>Recommended dilutions</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range: 80-160
<b>Specificity</b>	RAB30 Polyclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Ras-related protein Rab-30
<b>Gene Name</b>	RAB30
<b>Cellular localization</b>	Membrane ; Lipid-anchor ; Cytoplasmic side . Golgi apparatus, trans-Golgi network . Cytoplasm . Golgi apparatus .
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	22kD
<b>Human Gene ID</b>	27314
<b>Human Swiss-Prot Number</b>	Q15771
<b>Alternative Names</b>	
<b>Background</b>	similarity:Belongs to the small GTPase superfamily. Rab family.,



**ELK Biotechnology**



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night