



RHG28 rabbit pAb

Cat No.:ES10162

For research use only

Overview

Product Name	RHG28 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 250-330
Specificity	RHG28 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Rho GTPase-activating protein 28 (Rho-type GTPase-activating protein 28)
Gene Name	ARHGAP28 KIAA1314
Cellular localization	nucleoplasm,cytosol,integral component of membrane,cell junction,
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	80kD
Human Gene ID	
Human Swiss-Prot Number	Q9P2N2
Alternative Names	
Background	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1



ELK Biotechnology

Rho-GAP domain.,tissue specificity:Expressed in testis. Expressed at moderate level in kidney and ovary, and weakly expressed in spleen and skeletal muscle.,