

RL22L rabbit pAb

Cat No.:ES10181

For research use only

Overview

Product Name	RL22L rabbit pAb
Host species	Rabbit
Applications	WB;IP;ELISA;IP
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	IP 1:50-200 WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	RL22L Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	60S ribosomal protein L22-like 1
Gene Name	RPL22L1
Cellular localization	ribosome,cytosolic large ribosomal subunit,
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	13kD
Human Gene ID	200916
Human Swiss-Prot Number	Q6P5R6
Alternative Names	
Background	similarity:Belongs to the ribosomal protein L22e family.,