



ELK Biotechnology

STK38 Rabbit pAb
Catalog Number: ES10204
For research use only.

Overview

Product Name	STK38 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/5000-1/20000
Immunogen	Synthesized peptide derived from human protein . at AA range; 220-300
Specificity	STK38 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	STK38
Gene Name	STK38 NDR1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	51kDa
Human Gene ID	11329
Human Swiss-Prot Number	Q15208
Alternative Names	