



NEK10 rabbit pAb

Cat No.:ES10208

For research use only

Overview

Product Name	NEK10 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 270-350
Specificity	NEK10 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Serine/threonine-protein kinase Nek10 (EC 2.7.11.1) (Never in mitosis A-related kinase 10) (NimA-related protein kinase 10)
Gene Name	NEK10
Cellular localization	nucleus,cytoplasm,centrosome,protein kinase complex,
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	128kD
Human Gene ID	152110
Human Swiss-Prot Number	Q6ZWH5
Alternative Names	
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,cofactor:Magnesium.,similarity:Belongs to the protein kinase superfamily. NEK Ser/Thr protein kinase family. NIMA subfamily.,similarity:Contains 1 ARM repeat.,similarity:Contains 1 protein kinase domain.,

