



S6A13 rabbit pAb

Cat No.:ES10271

For research use only

Overview

Product Name	S6A13 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	S6A13 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Sodium- and chloride-dependent GABA transporter 2 (GAT-2) (Solute carrier family 6 member 13)
Gene Name	SLC6A13 GAT2
Cellular localization	Cell membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	66kD
Human Gene ID	6540
Human Swiss-Prot Number	Q9NSD5
Alternative Names	
Background	function:Terminates the action of GABA by its high affinity sodium-dependent reuptake into presynaptic terminals.,similarity:Belongs to the sodium:neurotransmitter symporter (SNF) family.,tissue specificity:Expressed in brain, kidney, lung, liver and testis.,

