



NAGT1 rabbit pAb

Cat No.:ES10292

For research use only

Overview

Product Name	NAGT1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 400-480
Specificity	NAGT1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Sodium-dependent glucose transporter 1
Gene Name	NAGLT1 KIAA1919 HSPC100
Cellular localization	Apical cell membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	56kD
Human Gene ID	91749
Human Swiss-Prot Number	Q5TF39
Alternative Names	
Background	function:May function as a sodium-dependent glucose transporter.,similarity:Belongs to the major facilitator superfamily.,

