



SUGP1 rabbit pAb

Cat No.:ES10307

For research use only

Overview

Product Name	SUGP1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SUGP1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	SURP and G-patch domain-containing protein 1 (RNA-binding protein RBP) (Splicing factor 4)
Gene Name	SUGP1 SF4
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	70kD
Human Gene ID	57794
Human Swiss-Prot Number	Q8IWZ8
Alternative Names	
Background	SF4 is a member of the SURP family of splicing factors.[supplied by OMIM, Sep 2003],

