



TCAL7 rabbit pAb

Cat No.:ES10369

For research use only

Overview

Product Name	TCAL7 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein AA range: 1-50
Specificity	TCAL7 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Transcription elongation factor A protein-like 7 (TCEA-like protein 7) (Transcription elongation factor S-II protein-like 7)
Gene Name	TCEAL7
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	11kD
Human Gene ID	56849
Human Swiss-Prot Number	Q9BRU2
Alternative Names	
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the TFS-II family. TFA subfamily.,

