

Cleaved-Lamin A (D230) rabbit pAb

Cat No.: ES1041

For research use only

Overview

Product Name Cleaved-Lamin A (D230) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human Lamin A. AA

range:181-230

Specificity Cleaved-Lamin A (D230) Polyclonal Antibody detects

endogenous levels of fragment of activated Lamin A protein resulting from cleavage adjacent to D230.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Prelamin-A/C

Gene Name LMNA

Cellular localization Nucleus . Nucleus envelope . Nucleus lamina.

Nucleus, nucleoplasm. Nucleus matrix.

Farnesylation of prelamin-A/C facilitates nuclear envelope targeting and subsequent cleavage by ZMPSTE24/FACE1 to remove the farnesyl group

produces mature lamin-A/C, whic

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 28+75kD
Human Gene ID 4000
Human Swiss-Prot Number P02545

+86-27-59760950

Alternative Names LMNA; LMN1; Prelamin-A/C



ELKbio@ELKbiotech.com

www.elkbiotech.con




Background

lamin A/C(LMNA) Homo sapiens The nuclear lamina consists of a two-dimensional matrix of proteins located next to the inner nuclear membrane. The lamin family of proteins make up the matrix and are highly conserved in evolution. During mitosis, the lamina matrix is reversibly disassembled as the lamin proteins are phosphorylated. Lamin proteins are thought to be involved in nuclear stability, chromatin structure and gene expression. Vertebrate lamins consist of two types, A and B. Alternative splicing results in multiple transcript variants. Mutations in this gene lead to several diseases: Emery-Dreifuss muscular dystrophy, familial partial lipodystrophy, limb girdle muscular dystrophy, dilated cardiomyopathy, Charcot-Marie-Tooth disease, and Hutchinson-Gilford progeria syndrome. [provided by RefSeq, Apr 2012],

HeLa
(kD)
11785483426-


Western Blot analysis of various cells using Cleaved-Lamin A (D230) Polyclonal Antibody



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using Lamin A (Cleaved-Asp230) Antibody. The picture on the right is blocked with the synthesized peptide.







(kD)

Western blot analysis of lysates from HeLa cells, treated with Etoposide 25uM 60', using Lamin A (Cleaved-Asp230) Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950