



RPRD2 rabbit pAb

Cat No.:ES10450

For research use only

Overview

Product Name	RPRD2 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 1370-1450
Specificity	RPRD2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Regulation of nuclear pre-mRNA domain-containing protein 2
Gene Name	RPRD2 KIAA0460 HSPC099
Cellular localization	nucleoplasm,DNA-directed RNA polymerase II, holoenzyme,
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	160kD
Human Gene ID	23248
Human Swiss-Prot Number	Q5VT52
Alternative Names	
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 CID domain.,

