

MAGD4 rabbit pAb

Cat No.:ES10556

For research use only

Overview

Product Name	MAGD4 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 300-380
Specificity	MAGD4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Melanoma-associated antigen D4 (MAGE-D4 antigen) (MAGE-E1 antigen)
Gene Name	MAGED4 KIAA1859 MAGED4A MAGEE1; MAGED4B
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	81kD
Human Gene ID	728239
Human Swiss-Prot Number	Q96JG8
Alternative Names	
Background	MAGE family member D4B(MAGED4B) Homo sapiens This gene is a member of the MAGED gene family. It is expressed only in brain and ovary, and some transcript variants of this gene are specifically expressed in glioma cells. This gene is clustered with other MAGED genes on chromosome Xp11. Alternatively spliced transcript variants encoding different isoforms have been described for





ELK Biotechnology

this gene. [provided by RefSeq, May 2011],



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C