



# ELK Biotechnology

**MOB1A Rabbit pAb**  
**Catalog Number: ES10654**  
**For research use only.**

## Overview

---

<b>Product Name</b>	MOB1A polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/500-1/2000</b> <b>ELISA</b> <b>1/5000-1/20000</b>
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range; 160-240
<b>Specificity</b>	MOB1A Polyclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	MOB1A
<b>Gene Name</b>	MOB1A C2orf6 MOB4B MOBK1B MOBKL1B
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>23kDa</b>
<b>Human Gene ID</b>	55233
<b>Human Swiss-Prot Number</b>	Q9H8S9
<b>Alternative Names</b>	