

EMID2 rabbit pAb

Cat No.:ES10685

For research use only

Overview

Product Name	EMID2 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 120-200
Specificity	EMID2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Collagen alpha-1(XXVI) chain (EMI domain-containing protein 2) (Emilin and multimerin domain-containing protein 2) (Emu2)
Gene Name	EMID2 COL26A1 EMU2
Cellular localization	Secreted, extracellular space, extracellular matrix.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	48kD
Human Gene ID	136227
Human Swiss-Prot Number	Q96A83
Alternative Names	
Background	This gene encodes a protein containing an emilin domain and two collagen stretches. This gene may be associated with aspirin-intolerant asthma. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2013],

