

Cleaved-MMP-1 22k (F100) rabbit pAb

Cat No.:ES1069

For research use only

Overview

Product Name	Cleaved-MMP-1 22k (F100) rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not
	yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human MMP1. AA
	range:81-130
Specificity	Cleaved-MMP-1 22k (F100) Polyclonal Antibody
	detects endogenous levels of fragment of activated
	MMP-1 22k protein resulting from cleavage adjacent
	to F100.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Interstitial collagenase
Gene Name	MMP1
Cellular localization	Secreted, extracellular space, extracellular matrix .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	19kD
Human Gene ID	4312
Human Swiss-Prot Number	P03956
Alternative Names	MMP1; CLG; Interstitial collagenase; Fibroblast
	collagenase; Matrix metalloproteinase-1; MMP-1
Background	matrix metallopeptidase 1(MMP1) Homo sapiens
	This gene encodes a member of the peptidase M10
	family of matrix metalloproteinases (MMPs).
	Proteins in this family are involved in the breakdown



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. The encoded preproprotein is proteolytically processed to generate the mature protease. This secreted protease breaks down the interstitial collagens, including types I, II, and III. The gene is part of a cluster of MMP genes on chromosome 11. Mutations in this gene are associated with chronic obstructive pulmonary disease (COPD). Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed. [provided by RefSeq, Jan 2016],

Western Blot analysis of various cells using Cleaved-MMP-1 22k (F100) Polyclonal Antibody



A549

(kD)

117-85-

Western blot analysis of lysates from A549 cells, treated
with etoposide 25uM 24h, using MMP1 (Cleaved-Phe100)
Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C