

IFNA4 rabbit pAb

Cat No.:ES10706

For research use only

Overview

Product Name	IFNA4 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of
	human protein AA range: 87-137
Specificity	IFNA4 Polyclonal Antibody detects endogenous
	levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Interferon alpha-4 (IFN-alpha-4) (Interferon
	alpha-4B) (Interferon alpha-76) (Interferon
	alpha-M1)
Gene Name	IFNA4
Cellular localization	Secreted.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	20kD
Human Gene ID	3441
Human Swiss-Prot Number	P05014
Alternative Names	
Background	function:Produced by macrophages, IFN-alpha have
	antiviral activities. Interferon stimulates the
	production of two enzymes: a protein kinase and an
	oligoadenylate synthetase.,polymorphism:Two
	forms exist; alpha-4a (shown here) and alpha-4b.
	They seem to be equally
	abundant., similarity: Belongs to the alpha/beta




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





interferon family.,

Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C