



TBD rabbit pAb

Cat No.:ES10739

For research use only

Overview

Product Name	TBD rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TBD Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Tubulin delta chain (Delta-tubulin)
Gene Name	TUBD1 TUBD
Cellular localization	Nucleus . Cytoplasm . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole . Cell projection, cilium . Associated with centrioles. Both cytoplasmic and nuclear. In the elongating spermatid it is associated with the manchette, a specialized microtubule system present during reshaping of the sperm head. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	49kD
Human Gene ID	51174
Human Swiss-Prot Number	Q9UJT1
Alternative Names	
Background	alternative products:Experimental confirmation may be lacking for some isoforms,function:In the elongating spermatid it is associated with the





ELK Biotechnology

manchette, a specialized microtubule system present during reshaping of the sperm head.,similarity:Belongs to the tubulin family.,subcellular location:Associated with centrioles. Both cytoplasmic and nuclear.,



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C