



# ELK Biotechnology

## Caspase-6 Rabbit pAb Catalog Number: ES1074 For research use only.

### Overview

---

Product Name	Caspase-6 polyclonal antibody
Host species	Rabbit
Applications	<b>WB,ELISA</b>
Species Cross-Reactivity	<b>Human</b>
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>ELISA:1/10000</b>
Immunogen	Synthesized peptide derived from the N-terminal region of human Caspase-6. at AA rangle; 20-100
Specificity	Caspase-6 Polyclonal Antibody detects endogenous levels of Caspase-6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Caspase-6
Gene Name	CASP6
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>33kDa</b>
Human Gene ID	839
Human Swiss-Prot Number	P55212
Alternative Names	CASP6; MCH2; Caspase-6; CASP-6; Apoptotic protease Mch-2