



ELK Biotechnology

Myc-tag Rabbit pAb
Catalog Number: ES1075
For research use only.

Overview

Product Name	Myc-tag polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	N/A
Recommended dilutions	Western Blot:1/1000-1/3000 ELISA:1/20000
Immunogen	Myc synthetic peptide conjugated to KLH.
Specificity	Myc-tag Polyclonal Antibody detects Myc-tagged recombinant proteins or Myc-tagged proteins overexpressed in cells.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.05% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Myc-tag
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	35kDa
Human Gene ID	
Human Swiss-Prot Number	
Alternative Names	c Myc; c Myc Epitope Tag; c Myc Tag; Myc Epitope Tag; Myc proto-oncogene protein; Transcription factor p64