



ELK Biotechnology

H31 Rabbit pAb
Catalog Number: ES10841
For research use only.

Overview

| | |
|---------------------------------|--|
| Product Name | H31 polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,ELISA |
| Species Cross-Reactivity | Human,Mouse,Rat |
| Recommended dilutions | Western Blot:1/500-1/2000 ELISA 1/5000-1/20000 |
| Immunogen | |
| Specificity | H31 Polyclonal Antibody detects endogenous levels of protein. |
| Formulation | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | H31 |
| Gene Name | HIST1H3A H3FA; HIST1H3B H3FL; HIST1H3C H3FC; HIST1H3D H3FB; HIST1H3E H3FD; HIST1H3F H3FI; HIST1H3G H3FH; HIST1H3H H3FK; HIST1H3I H3FF; HIST1H3J H3FJ |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 14kDa |
| Human Gene ID | 8350 |
| Human Swiss-Prot Number | P68431 |
| Alternative Names | |