



OLIG3 rabbit pAb

Cat No.:ES10942

For research use only


Overview

Product Name	OLIG3 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	OLIG3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Oligodendrocyte transcription factor 3 (Oligo3) (Class B basic helix-loop-helix protein 7) (bHLHb7) (Class E basic helix-loop-helix protein 20) (bHLHe20)
Gene Name	OLIG3 BHLHB7 BHLHE20
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	29kD
Human Gene ID	167826
Human Swiss-Prot Number	Q7RTU3
Alternative Names	
Background	function:May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress specification of class B neurons.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,





ELK Biotechnology



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.