

## SG3A2 rabbit pAb

Cat No.: ES11168

For research use only

## Overview

Product Name SG3A2 rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein .

at AA range: 20-100

**Specificity** SG3A2 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

StorageStore at -20°C. Avoid repeated freeze-thaw cycles.Protein NameSecretoglobin family 3A member 2 (Pneumo

secretory protein 1) (PnSP-1) (Uteroglobin-related

protein 1)

Gene Name SCGB3A2 PNSP1 UGRP1 UNQ566/PRO1128

Cellular localization Secreted.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 10kD
Human Gene ID 117156
Human Swiss-Prot Number Q96PL1

**Alternative Names** 

**Background** The protein encoded by this gene is a secreted lung

surfactant protein and a downstream target of thyroid transcription factor. A single nucleotide polymorphism in the promoter of this gene results in susceptibility to asthma.[provided by RefSeq, Mar


2010],



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com





Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night

