



OTUD4 rabbit pAb

Cat No.:ES11195

For research use only

Overview

Product Name	OTUD4 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 910-990
Specificity	OTUD4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	OTU domain-containing protein 4 (HIV-1-induced protein HIN-1)
Gene Name	OTUD4 HIN1 KIAA1046
Cellular localization	Cytoplasm . Nucleus . Primarily cytoplasmic. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	122kD
Human Gene ID	54726
Human Swiss-Prot Number	Q01804
Alternative Names	
Background	Alternatively spliced transcript variants have been found for this gene. The smaller protein isoform encoded by the shorter transcript variant is found only in HIV-1 infected cells. [provided by RefSeq, Jul 2010],

