



ERR2 rabbit pAb

Cat No.:ES11197

For research use only

Overview

Product Name	ERR2 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 270-350
Specificity	ERR2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Steroid hormone receptor ERR2 (ERR beta-2) (Estrogen receptor-like 2) (Estrogen-related receptor beta) (ERR-beta) (Nuclear receptor subfamily 3 group B member 2)
Gene Name	ESRRB ERB2 ESRL2 NR3B2
Cellular localization	Nucleus . Cytoplasm . Chromosome .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	55kD
Human Gene ID	2103
Human Swiss-Prot Number	O95718
Alternative Names	
Background	This gene encodes a protein with similarity to the estrogen receptor. Its function is unknown; however, a similar protein in mouse plays an essential role in placental development. [provided by RefSeq, Jul 2008],

