

## **EPYC** rabbit pAb

## Cat No.:ES11217

For research use only

## Overview

Product Name	EPYC rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein .
	at AA range: 140-220
Specificity	EPYC Polyclonal Antibody detects endogenous levels
	of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Epiphycan (Dermatan sulfate proteoglycan 3)
	(Proteoglycan-Lb) (PG-Lb) (Small
	chondroitin/dermatan sulfate proteoglycan)
Gene Name	EPYC DSPG3 PGLB SLRR3B
<b>Cellular localization</b>	Secreted, extracellular space, extracellular matrix .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	35kD
Human Gene ID	1833
Human Swiss-Prot Number	Q99645
Alternative Names	
Background	epiphycan(EPYC) Homo sapiens Dermatan
	sulfate proteoglycan 3 is a member of the small
	leucine-rich repeat proteoglycan family. This gene is
	composed of seven exons. It regulates fibrillogenesis
	by interacting with collagen fibrils and other
	extracellular matrix proteins. [provided by RefSeq,
	Jul 2008],




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C