



# ELK Biotechnology

**KLK15 Rabbit pAb**  
**Catalog Number: ES11229**  
**For research use only.**

## Overview

---

Product Name	KLK15 polyclonal antibody
Host species	Rabbit
Applications	<b>WB,ELISA</b>
Species Cross-Reactivity	<b>Human</b>
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>ELISA</b> <b>1/5000-1/20000</b>
Immunogen	Synthesized peptide derived from human protein . at AA range; 20-100
Specificity	KLK15 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	KLK15
Gene Name	KLK15
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>28kDa</b>
Human Gene ID	55554
Human Swiss-Prot Number	Q9H2R5
Alternative Names	